


PROTAPER • GOLD™



The same winning technique with **greater flexibility**



THEN: PROTAPER® UNIVERSAL sets a new standard in efficiency



NOW: The PROTAPER® legacy continues with PROTAPER GOLD™

SAME: File Sequence • File Sizes • Motor Settings

DIFFERENT: Heat Treatment • Gold Colour

IMPROVED: Flexibility • Cyclic Fatigue Resistance • Tactile Feeling • Safety

**+
WE
KNOW
ENDO.**

**PROTAPER GOLD™
IS NOW REPLACING
PROTAPER® UNIVERSAL**

PROTAPER•GOLD™

Product references

PROTAPER GOLD™ FILES



SX STERILE
19MM A 0410 219 G01 03



S1 STERILE
21MM A 0410 221 G01 03
25MM A 0410 225 G01 03
31MM A 0410 231 G01 03



F3 STERILE
21MM A 0411 221 G03 03
25MM A 0411 225 G03 03
31MM A 0411 231 G03 03



S2 STERILE
21MM A 0410 221 G02 03
25MM A 0410 225 G02 03
31MM A 0410 231 G02 03



F4 STERILE
21MM A 0411 221 G04 03
25MM A 0411 225 G04 03
31MM A 0411 231 G04 03



F1 STERILE
21MM A 0411 221 G01 03
25MM A 0411 225 G01 03
31MM A 0411 231 G01 03



F5 STERILE
21MM A 0411 221 G05 03
25MM A 0411 225 G05 03
31MM A 0411 231 G05 03



F2 STERILE
21MM A 0411 221 G02 03
25MM A 0411 225 G02 03
31MM A 0411 231 G02 03

ASSORTMENT FILE PACKS
21MM SX/F3 A 0409 221 G01 03
25MM SX/F3 A 0409 225 G01 03
31MM SX/F3 A 0409 231 G01 03



“PROTAPER GOLD™ produces consistently predictable shapes for predictable cleaning and obturation that is easier and safer than ever before.”

Dr. John West²

All current PROTAPER® Universal obturation methods can be used with PROTAPER GOLD™.

In some instances PROTAPER GOLD™ files may appear slightly curved. This feature is due to the GOLD technology allowing the files to be pre bent to achieve better access.

PROTAPER GOLD™ now features an 11mm handle for improved intra oral visibility (versus 13mm with PROTAPER® Universal).



For more information, please visit www.dentsplymaillefer.com

1. Average value for PROTAPER GOLD™ compared to PROTAPER® Universal. Based on internal testing. Data on file.
2. PROTAPER GOLD™ (PTG) Shapers and Finishers. Dentsply Sirona (DTDS).